



MARKING

3-0231



Abdullayev Corneal Marker for Keratoplasty



- 10.00 & 11.00 mm diameters
- With central marking point
- Improves centration of the cornea
- Allows quick placement of the central dot

2-0331T

Grooved Mendez Degree



- Used to determine incision angles
- With 4 grooves
- 3.00/16.50 mm diameters

3-143T

Tunnel Zone Marker




- Marking tunnel for ring position
- Initial lines for incisions
- 4.00/6.00 mm diameter
- 2.00 mm lines

TUNNEL CREATION

6-00/6-020

Astigmatic and Arcuate Keratotomy Diamond Knife



- 1.00 mm radial incision at 70–80% of corneal thickness
- Micrometric scale to customize the depth of the blade from 0 to 1000 microns (0 to 1.00 mm)
- Increment – 5.00 microns (0.005 mm)

**16-175
16-176**

LWVI CAIRS Dissectors ✦



NEW

- Designed for accuracy, comfort, and versatility
- 7 sizes for both left and right approach
- Blade widths: 1.25 / 1.50 / 2.00 / 2.50 mm

Inner diameter 4.50 mm

16-175-1.25LS	16-175-1.25RS
16-175-1.50LS	16-175-1.50RS
16-175-2.00LS	16-175-2.00RS
16-175-2.50LS	16-175-2.50RS

Inner diameter 5.50 mm

16-176-1.25LS	16-176-1.25RS
16-176-1.50LS	16-176-1.50RS
16-176-2.00LS	16-176-2.00RS



16-315T

LWVI CAIRS Vacuum Stabilization Guide ✦

Suction Ring equipped with a vacuum connection port:

- Ensures stable centering on the cornea
- Keeps dissectors firmly positioned during intrastromal tunnel formation



10-035

Elevator for ICSR Implantation




- To create a small pocket in intended direction of the tunnel with a mild rotatory movement

TUNNEL CREATION

13-146


Suarez Spreader



- To enlarge the small corneal pockets in both directions in a proper depth

13-147

Bicalto Guide



- Allows easy Tunnel Maker insertion into tunnel

BASIC

16-036T

Fine Thornton Phaco Fixation Ring with Swivel




- 3/4 open 13.00 mm diameter ring with teeth

18-304

Plastic Sterilization Tray with Silicone Finger Mat, Extra Large

- Sterilizable by all standard methods




SEGMENTS IMPLANTATION



NEW

5-080
5-081

LWVI CAIRS Y-Shaped Pusher, Right ✦
LWVI CAIRS Y-Shaped Pusher, Left ✦



- Uniquely curved “Y” design
- Enhanced reach, control, and ease of segment placement

4-0822T

Castroviejo Fixation Forceps




0.50 mm, 1x2 teeth



- For holding the globe in place
- Handling and fixation of conjunctiva and episcleral tissue

4-178S


McPherson Straight Tying Forceps



- Hard-faced with tungsten carbide coating
- Greater grip and increased control
- For 8-0 and 11-0 sutures
- 7.00 mm tying platforms

5-0322

Reversed Sinskey Hook



- For easy stromal rings manipulation